



**PHILIPS**


HR2621, HR2650  
HR2651, HR2652  
HR2653, HR2655  
HR2656, HR2657



		MAX.	⌚	📊
		100-200 g	30 sec.	
		100-400 ml	60 sec.	
		100-500 ml	60 sec.	
		100-1000 ml	60 sec.	
		250 ml	70-90 sec.	
		4 x	120 sec.	
		250 ml	70-90 sec.	
		4 x	120 sec.	
		750 g	180 sec.	
		750 g	120 sec.	
		200 g	5 x 1 sec.	
		200 g	5 sec.	
		200 g	10 sec.	
		30 g	10 sec.	
		200 g	15 sec.	
		200 g	30 sec.	
		100 g	20 sec.	
		80 g	30 sec.	




HR2621, HR2652  
HR2655, HR2657



HR2653



HR2621, HR2652  
HR2653, HR2655, HR2657




**HR2656**  
**HR2657**




**HR2650, HR2651, HR2652**  
**HR2653, HR2655, HR2656, HR2657**



**HR2651**  
**HR2655**



**HR2651**  
**HR2655**



**HR2656**  
**HR2657**

	1	2	3	

	MAX.		
		400 ml	60 sec.
		400 ml	60 sec.

	MAX.		
		20 min	
		60 sec.	

**HR2621**  
**HR2652**  
**HR2655**  
**HR2657**

	✗	✗	✓
	✓	✗	✓
	✗	✗	✓
	✓	✓	✓
	✗	✗	✓
	✓	✓	✓
	✗	✗	✓
	✓	✓	✓
	✓	✓	✓
	✓	✓	✓
	✓	✓	✓
	✓	✓	✓

